

Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US05/005616

International filing date: 23 February 2005 (23.02.2005)

Document type: Certified copy of priority document

Document details: Country/Office: US
Number: 60/547,379
Filing date: 23 February 2004 (23.02.2004)

Date of receipt at the International Bureau: 23 March 2005 (23.03.2005)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse

1294646

THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

March 10, 2005

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.

APPLICATION NUMBER: 60/547,379

FILING DATE: February 23, 2004

RELATED PCT APPLICATION NUMBER: PCT/US05/05616



Certified by

Under Secretary of Commerce
for Intellectual Property
and Director of the United States
Patent and Trademark Office

17157 U.S. PTO
022304

PTO/SB/16 (01-04)
Approved for use through 07/31/2006. OMB 0651-0032
U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE
Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

PROVISIONAL APPLICATION FOR PATENT COVER SHEET

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

Express Mail Label No. **ET 590117968 US**

22582 U.S. PTO
60/547379

022304

INVENTOR(S)						
Given Name (first and middle [if any])		Family Name or Surname		Residence (City and either State or Foreign Country)		
EUDES FRANCOIS NARIE		de CRECY		LA JOLLA, CA, 92037		
Additional inventors are being named on the _____ separately numbered sheets attached hereto						
TITLE OF THE INVENTION (500 characters max)						
CONTINUOUS CULTURE APPARATUS WITH MOBILE VESSEL, ALLOWING SELECTION OF FETTER CELL VARIANTS						
Direct all correspondence to: CORRESPONDENCE ADDRESS						
<input type="checkbox"/> Customer Number: _____						
OR						
<input checked="" type="checkbox"/> Firm or Individual Name		EUDES FRANCOIS NARIE de CRECY				
Address		5919 LA JOLLA CORONA DRIVE				
Address						
City		LA JOLLA	State	CA	Zip	92037
Country		USA	Telephone	(858) 454 8758	Fax	(858) 454 8758
ENCLOSED APPLICATION PARTS (check all that apply)						
<input checked="" type="checkbox"/> Specification Number of Pages		2		<input type="checkbox"/> CD(s), Number _____		
<input checked="" type="checkbox"/> Drawing(s) Number of Sheets		2		<input type="checkbox"/> Other (specify) _____		
<input type="checkbox"/> Application Data Sheet. See 37 CFR 1.76						
METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT						
<input type="checkbox"/> Applicant claims small entity status. See 37 CFR 1.27.				FILING FEE Amount (\$)		
<input checked="" type="checkbox"/> A check or money order is enclosed to cover the filing fees.				80.00		
<input type="checkbox"/> The Director is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number: _____						
<input type="checkbox"/> Payment by credit card. Form PTO-2038 is attached.						
The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.						
<input checked="" type="checkbox"/> No.						
<input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are: _____						

Respectfully submitted,

[Page 1 of 2]

Date **02/23/04**

SIGNATURE **EUDES DE CRECY**

REGISTRATION NO. _____

TYPED or PRINTED NAME **EUDES DE CRECY**

(if appropriate)

TELEPHONE **(858) 454 8758**

Docket Number: _____

USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT

This collection of information is required by 37 CFR 1.51. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Mail Stop Provisional Application, Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

PROVISIONAL APPLICATION COVER SHEET
Additional Page




PTO/SB/16 (08-03)

Approved for use through 07/31/2008. OMB 0651-0032

U.S. Patent and Trademark Office, U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

Docket Number

INVENTOR(S)/APPLICANT(S)		
Given Name (first and middle [if any])	Family or Surname	Residence (City and either State or Foreign Country)
EUDES FRANÇOIS NARIE	de CRECY	LA JOLLA, CA, 92037
		
		
		

[Page 2 of 2]

Number _____ of _____

WARNING: Information on this form may become public. Credit card information should not be included on this form. Provide credit card information and authorization on PTO-2038.

Title: Continuous culture apparatus with mobile vessel, allowing selection of fitter cell variants.

SUMMARY

DESCRIPTION

Continuous culture vessels have been used since the 1950's to select for cells that have a higher proliferation rate under given conditions (Monod J., 1950, Ann. Inst. Pasteur. 19:390-410; Novick A., Szilard L., 1950, Science 112:715-716). The power of these devices have never been used to their fullest because they positively select mutants that evade dilution, primarily through attachment to vessel surfaces, resulting in persistent sub-populations of uncontrollable size and growth rate.

One possible solution to overcome this drawback is the implementation of a device with two growth chambers periodically undergoing transient phases of sterilization as described in the patent by Marliere and Mutzel (1999, DE2982162U1).

The solution described in this application is radically different as it consists of using a mobile vessel around the solution instead of moving the solution from one chamber to another.

More specifically, the vessel consists of a sterile flexible tube containing the culture media.

A portion of this tube that can be physically and temporarily separated from the rest of the tube by some lock devices, will contain the cell culture, and will be designated therein as the culture vessel. Fresh medium will be periodically added to the culture by sliding a portion of medium filled tubing.

The design of the lock, see figure, will drive the amount of medium that will be added at each sliding cycle.

Spend culture and medium will also be eliminated through a tubing sliding step.

The frequency of the movement of the tube can be set by the experimenter reproducing a chemostat regime or can be for example regulated through a feedback loop depending on the turbidity of the culture in the vessel that will be read via a turbidometer, emulating a turbidostat regime. Other regulation approaches can be implemented.

The system is designed to allow the complete rotation of the vessel and the locks.

The positioning of the vessel will be adjusted in accordance to the aeration cycle.

Aeration can be provided two ways:

- 1- by passive diffusion of ambient air through gas permeable tubing,
- 2- by injecting filtered air that will first bubble through the media and then be introduced in the vessel through the entrance locking sub-system.

Previously the spend air will have been eliminated by rotation of the vessel in the upright position as described in attached figures.

The presence of several gates (superior to 2) should avoid any contamination in the tubing of the medium in the up chamber tubing. However several UV gates can be added upstream and downstream of the culture vessel for additional security.

Inoculation will be performed by sterily injecting a cell culture through the lower lock that would have been set in the opened position.


An added lock will eventually be added downstream of the chamber lock to create a sampling chamber.

Tubing specifications: The tubing must be large enough to avoid capillarity problems that would prevent the air from circulating freely. The tubing should sustain sterilization cycle, be gas permeable, be flexible enough to be closed by compression and be clear enough to allow any optical measurement.

CONTINUOUS CULTURE APPARATUS WITH MOBILE VESSEL
ALLOWING SELECTION OF RITER CELL VARIANTS

DRAWING N°1

AIR FILTER



CONTROL
SYSTEM

LOCK 1

LOCK 2

LOCK 3







LOCK 4

LOCK 5

LOCK 6

LOCK 7

LOCK 8



E1 = ENGINE 1

E2 = ENGINE 2


T = TURBIDOMETER

UN = ULTRA-VIOLET GATE


V.D. = VOLUME OF DILUTION

LOCK TOOTH DETAIL

DRAWING NO. 2



POSITION 2



POSITION 1